


PikoDes P



PikoDes P, a powerful, highly concentrated oxidizing powder disinfectant to destroy fungi, viruses, bacteria and yeasts in an effective way.

Benefits:

- **Very fast working and active at low concentrations**
- **Extremely effective against specific pathogenic micro-organisms**
- **Based on multiple active substances with a very broad spectrum of activity**
- **PikoDes P is completely biodegradable**



Suitable for:

- Disinfection of hard surfaces and materials in places where food and drinks are prepared, handled and stored
- Disinfection of surfaces, materials and tools in animal housings, boots etc.

How to use:

- Clean the surface thoroughly by using a powerful soaking product (e.g. PikoClean Mousse LA Sense)
- Dosage 1-2%
- Apply PikoDes P by spraying* (make sure that no unprotected animals or persons are present while spraying)
- Make sure the surface remains wet during the exposure time
- After soaking (5-30 minutes) rinse thoroughly with water

Effective against:

Fungi and yeasts: Candida albicans, Microsporum canis, Trichophyton mentagrophytes

Bacteria: Bacillus cereus, Campylobacter jejuni, Clostridium perfringens, Dermatophilus congolensis, Enterococcus faecalis, Enterococcus hirae, Escherichia coli, Haemophilus equigenitalium, Helicobacter pylori, Listeria monocytogenes, Proteus vulgaris, Pseudomonas aeruginosa, Salmonella choleraesuis, Salmonella typhimurium, Staphylococcus aureus, Streptococcus equi, Vibrio cholerae

Viruses: Aujeszky's disease, Avian Influenza H3N2, Bovine Pneumonia, Chicken anemia, Circo virus, Egg Drop Syndrome, Enteritis, Foot and Mouth disease, Gumboro, Infectious Bovine Rhinotracheitis, Infectious Pancreatic Necrosis, Infectious Salmon Anaemia, Newcastle, Parvo Viral Enteritis, Porcine Ear Blue/PRRS, Respiratory disease complex, Swine Vesicular disease, Tenosynovitis, Talfan Teschen disease, Viral diarrhoea, African Swine Fever (ASF), Enteric cytopathic bovine orphan virus (ECBO)

Contains:

Pentapotassium bis(peroxymonosulphate) bis(sulphate) 45%
Sodium dichloroisocyanurate 2,5%

Dosage:

1-2% (100g-200g at 10L water)

Soaking time:

5 - 30 minutes
(For viruses at least 30 minutes)

Storage:

Cool, dry

Pack size (kg):

10 | 1000

*PikoDes P is a biocide. Use PikoDes P (and all other biocides) safely. Always read the label and product information before use